



# TECHTRAINER™

Simply Intelligent.





# TECHTRAINER™



## A new standard in affordable, functional metrics

BTE Technologies developed the TechTrainer to provide physical therapists, athletic trainers, and strength conditioning specialists with a compact, powerful system that provides objective data capture with the best in pneumatic resistance technology. Never before have such advanced, easy-to-use measurements been integrated into a high performance pneumatic cable system.






## Motion

- Detect discrete differences
- Safe, smooth at any speed
- Baseline, and progress




### Upper Extremity



### Lower Extremity



### Core Conditioning




***“For me, the TechTrainer is marketable to the masses...  
We can evaluate an athlete’s strengths with a level of data quality that  
hasn’t been available before for smaller businesses and private practices...”***

—Tristan Tillette, CSCS  
Owner & Head Performance Coach  
Elite Sports Performance Birmingham, AL



# Intelligence

- Real time feedback
- Objective reports
- Client Motivation



Peak Power Test – Bilateral Comparison



Power Training (Sets)



# Powerful, Easy to Use Software

TechTrainer provides real time feedback for client motivation, pre-defined progressive evaluation & rehab programs, and comprehensive reports. The intuitive touchscreen software is easy to use and navigate for clinicians as well as clients.



Endurance Test – Bilateral Comparison



Personalized Testing and Training Options

*“Power is not just for athletes... It is a crucial component in any rehab program focused on bridging the gap between injury and high level function. This extremely powerful rehab and training tool should be in every clinic in the country. BTE has created one of the most exciting training and rehab products I have ever seen!”*

—Darcy Norman PT, ATC, CSCS  
Alpine, Endurance, and Performance Therapy Consultant  
Fitness and Rehab Coach,  
FC Bayern Munich

# Detailed printed reports



- Detailed Evaluations
- Summary Evaluations
- Treatment Progress

Capture incremental revenue through physical performance testing with written reports.

# Motion, Intelligence, Simplicity.

# TECHTRAINER™

BTE Technologies is the leader in testing and rehabilitation solutions. Our intuitive TechTrainer software reflects over 25 years of experience in designing functional technological solutions for healthcare and rehabilitation industries. Our systems provide objective data – so you can assess patients with greater precision, track progress with higher accuracy, and provide the type of outcomes that your referral sources and clients need. With over 5,000 installations worldwide, our technology is used by a broad variety of clinics and industries to improve their professional outcomes and quality of care. Call us today or visit us at [www.btetech.com](http://www.btetech.com).

## THE TRUE LEADER IN:

Physical  
Evaluation &  
Rehab

Employee  
Functional  
Capacity

Human  
Performance  
Assessment

Post-Offer of  
Employment  
Testing

Medical  
Impairment  
Assessment

## Simplicity

- Practical
- Functional
- Proven

- A** Powerful touchscreen computer keeps clients motivated with realtime data feedback, printable reports, and customizable testing/training progress tracking.
- B** High speed bilateral cable system with rotating pulleys for unrestricted natural movement through complete range of motion.
- C** Straps and Attachments to simulate numerous functional activities, including baseball, tennis conveniently store on durable, lightweight tool rack.
- D** Arms independently adjustable to 12 positions each, allow for replicating a multitude of functional motions.
- E** On-the-fly easy force adjustment in 1 lb increments allows safe use with clients from all walks of life.
- F** Quiet, clean air compressor powers up to six TechTrainers.
- G** Compact footprint as shown, or legs can be removed and system can be bolted to concrete floor for even more space savings.



web

[www.btetech.com](http://www.btetech.com)

phone

1-800-331-8845 [ U.S. and Canada ]  
1-410-850-0333 [ other countries ]



Corporate Headquarters

7455-L New Ridge Road  
Hanover, MD 21076

Experience. Technology.

© 2008 BTE Technologies, Inc. All rights reserved.