

**R001**

0.29 mJ/mm<sup>2</sup>

Penetration Depth (mm) per Energy Level	
Minimum (E1)	4
Maximum (E6)	21

**R005**

0.38 mJ/mm<sup>2</sup>

Penetration Depth (mm) per Energy Level	
Minimum (E1)	50
Maximum (E6)	59

**R020**

0.38 mJ/mm<sup>2</sup>

Penetration Depth (mm) per Energy Level	
Minimum (E1)	66
Maximum (E6)	75

**R035**

0.38 mJ/mm<sup>2</sup>

Penetration Depth (mm) per Energy Level	
Minimum (E1)	82
Maximum (E6)	91