



Nuovo bracciolo
brevetto Surace
New armrest design
Surace patented

➔ Standard ma robusta

|| Squillo Robusta sedia pieghevole

Standard: Telaio stabilizzante. Braccioli desk estraibili. Pedane girevoli ed estraibili. Freni.

Optional: Appoggiagambe regolabili, girevoli ed estraibili. Ruote ad estrazione rapida.

🇬🇧 Squillo Solid folding chair


Standard: Stabilised frame. Removable armrests and desk. Swinging, removable footrests. Brakes.

Optional: Adjustable, swinging and removable legrests. Quick release wheels.

Classificazione ISO 12.21.06.039



=



cm Ø 60	cm Ø 20	L 43	L 61	L 75	L 105	h 51	Standard	Optional	Portata 120 kg